




Insight into images

Robert Martí

robert.marti@udg.edu



- Insight into images
- Contents
 - Image information
 - Image Modalities
 - Image formats
 - DICOM, PACs

Some slides taken from O. Díaz (UB).

Before starting ...

What is a digital hospital?

Before starting ...

What sort of (**digital**) health data can you think of?

Introduction

- **Health Record** (non-imaging data)
 - Clinical document with patient information to provide support to healthcare members.
 - Includes all medical information regarding a person.
 - Unique ID number to each patient/client.
 - Confidential.
 - Potential problems:
 - Large pieces of information -> time.
 - Fragmentation.
 - Illegibility.
 - Anonymisation


} Handwriting

Welcome to the digital data era:
Electronic Health Record (EHR)


Introduction

- **Medical image** (imaging data)

Analogic




Digital




Introduction

- **Medical image (imaging data)**





2D



2.2D



3D




4D

Image Information


Pixels and images

- What is an image?
- Discrete grid representation of a continuous signal




Pixels and images

- Sampling grid



Pixels and images

- Image index



Pixels and images

- Pixel and physical coordinates
- Conversion

$$P[0] = \text{Index}[0] \times \text{Spacing}[0] + \text{Origin}[0]$$


$$P[1] = \text{Index}[1] \times \text{Spacing}[1] + \text{Origin}[1]$$

$$\text{Index}[0] = \text{floor} \left(\left(P[0] - \text{Origin}[0] \right) / \text{Spacing}[0] + 0.5 \right)$$

$$\text{Index}[1] = \text{floor} \left(\left(P[1] - \text{Origin}[1] \right) / \text{Spacing}[1] + 0.5 \right)$$


Pixels and images

- Image Region




Pixels and images

- Resampling: sub-sampling




Pixels and images

- Resampling: super-sampling



Pixels and images

- How do we compare images?




Pixels and images

- Comparing images




Patient A: 600 pixels lesion




Patient B: 800 pixels lesion


Pixels and Images




Patient A: 300mm² lesion



Pixels and Images



Patient A: lesió de 300mm^2







Pixels and Images





Image Modalities

Image Modalities




Breast imaging

- Mammography




Breast Imaging

- Compression
 - Minimize patient movement
 - Increase imaging area (minimise tissue overlap)
 - Less noise (reduce breast thickness, less scattering)





Breast imaging

- FFDM
 - 50-100 μm in pixel size
 - 12-16 bits per pixel
 - Image size of 8-50 MB




Breast imaging

- Phantoms




Breast imaging

- DBT Tomosynthesis




Breast Imaging

- Dynamic Contrast Enhanced MRI (DCE-MRI)




Breast imaging

- 3D ultrasound



Information in the hospital


- HIS. Hospital Information System.
 - Clinical patient information
 - Prescriptions
 - Operating room management
 - Patient management
 - Logistics
 - Financing
 - Public administration
 - Softwares: SAP, Costaisa, HP
- Protocol: HL7, DICOM, XML, Webservices, etc.
- RIS / PACS / DICOM ?



- Interoperability standard
- Technology and platform independent
- <http://www.hl7.org/>
- Patient information interchange standard based on special characters separation.

```
MSH|^~\&|EPIC|SYS|HOSP|ADT|201502031126|SEC|ADT^A01|001199|P|2.3
EVN|A01|201502031126
PID|||12001||SIMPSON^HOMER||19670824|M|||742 Evergreen Terrace St. ^^ Springfield ^ OR
^ 90020 ^ USA|||||
NK1|1|SIMPSON^MARGE|WIFE|||||NK
PV1|1||2000^2012^01||||11277^HIBBERT^JULIUS^J|||SUR||-||ADM|A0-
AL1|1||^Penicillin||Hives
```

- **PID**: patient ID
- **NK1**: next of kin
- **PV1**: patient visit
- **AL1**: allergies




- Integrating the Healthcare Enterprise
- Coordinated by the Radiological Society of North America (RSNA).
- Consortium between companies and public health organisations for improving interconnectivity in healthcare systems

<https://www.ihe.net/>




RIS / PACS / DICOM

- Radiology Information System (RIS)
- Picture Archiving and communication system (PAC)




DICOM

- DICOM: Digital Imaging and COmmunications in Medicine
- PACS: Picture Archiving and Communication System



DICOM & PACS



DICOM

- Digital Imaging and Communications in Medicine.
- Proposed by the American College of Radiology (ACR) and the National Electrical Manufacturers Association (NEMA)



<https://www.dicomstandard.org/>

























































































































- Standard to manage, store, print, and communicate medical imaging information and related data.
- Image format but also communication protocol

<http://www.dclunie.com/dicom-status/status.html>

DICOM origins

- 1982: ACR-NEMA Committee formed
- 1983: ACR-NEMA meeting
- 1985: ACR-NEMA version 1.0
 - 1st data saving in media and communication between non-proprietary soft
- 1988: ACR-NEMA version 2.0
 - Terminology, data structure and encoding.
- **1992: ACR-NEMA version 3.0 or DICOM**
 - Rely on ISO/OSI model and use of TCP/IP protocol.
 - Unique identifiers (GGGG, NNNN).
- Supplements
 - **Suppl. 28 Grayscale Standard Display Function (GSDF), 1998**
 - Zoom, windowing, or annotations

DICOM origins



Title	Format (see Key below)					
	PDF	HTML	CHTML	DOCX	ODT	XML
DICOM Part 1: Introduction and Overview						
DICOM Part 2: Conformance						
DICOM Part 3: Information Object Definitions						
DICOM Part 4: Service Class Specifications						
DICOM Part 5: Data Structures and Encoding						
DICOM Part 6: Data Dictionary						
DICOM Part 7: Message Exchange						
DICOM Part 8: Network Communication Support for Message Exchange						
DICOM Part 10: Media Storage and File Format for Media Interchange						
DICOM Part 11: Media Storage Application Profiles						
DICOM Part 12: Media Formats and Physical Media for Media Interchange						
DICOM Part 14: Grayscale Standard Display Function						
DICOM Part 15: Security and System Management Profiles						
DICOM Part 16: Content Mapping Resource						
DICOM Part 17: Explanatory Information						
DICOM Part 18: Web Services						
DICOM Part 19: Application Hosting						
DICOM Part 20: Imaging Reports using HL7 Clinical Document Architecture						
DICOM Part 21: Transformations between DICOM and other Representations						
DICOM Part 22: Real-Time Communication						
DICOM Parts 1-22: Bulk Download	Zip file for each format					

What's in the DICOM?

- Patient information
 - Patient: name, age , ...
 - Image acquisition: system, parameter, organ, who, where,
 - Image information: resolution, spacing, ...
 - Image or Images
 - ...
- But it is also a **communication** protocol between machines, PACs, workstations, etc.

What's in the DICOM?

- Hierarchical structure




DICOM file

File extension: .dcm


```
{0010,0010}: Anonymised00034  
{0010,0020}: XXX_72947  
{0010,0030}: 19000101  
{0010,0040}: M  
{0010,1010}: 039Y  
{0018,0015}: HEAD
```

Header



Image

DICOM file



Data Tag Element	16-bits unsigned representing Group and Element number.
Value Representation (VR)	2 bytes characters containing the VR of Data Element (e.g. string)
Value Length	<ul style="list-style-type: none"> - 16 or 32-bit (explicit or implicit) unsigned int containing length of the Value Field. - 32-bit length. Used for certain VR (SQ, UN, OW, OB)
Value Field	Even number of bytes containing the value(s) of Data Element

DICOM file

VR	Definition	Length	Sample
AS	Age string	4 bytes fixed	047Y
CS	Code string	16 bytes max.	MG
DA	Date string	8 bytes fixed	19930822
DS	Decimal string	16 bytes max.	44.8
LO	Long string	64 chars. max.	24X29
PN	Person Name	64 chars max	Simpson^Homer
SH	Short String	16 chars max.	H0711151700038
TM	Time	14 bytes max.	151850
UI	Unique Identifier	64 bytes max.	1.2.840.10008.1.2.1
US	Unsigned Short	2 bytes fixed	16

DICOM file

Tag / ID	Keyword	VR	Sample
(0008,0020)	StudyDate	DA	20160711
(0008,0022)	AcquisitionDate	DA	20160711
(0008,0060)	Modality	CS	MG
(0010,0010)	PatientName	PN	XXXXX_J9EKL8218SO7CHGZ
(0010,0020)	PatientID	LO	XXXXX_HX3LCLOW32528JVC
(0019,1025)	<HOLOGIC, Inc.> [25]	SH	FAST
(0028,0010)	Rows	US	4096
(0028,011)	Columns	US	3328
(0028,0030)	PixelSpacing	DS	0.065238 / 0.065238
(0028,0100)	BitsAllocated	US	16

DICOM file

- 8 bits or 16 bits?

```
1001101001110100110110011000011010110011010101100
0011001111001011001010101111000101000110011100010
1100101010011100101111000100011101010100001101010
1001001110011010011101001101100110011000010111011
1001010101110011010110011010101100001100111100101
1001010100100011011101011110001010001100111000101
1001010100111001011110001000111010101000011010101
```


...

$$2^8 = 256 \rightarrow [0-255]$$


$$2^{16} = 65,536 \rightarrow [0-65,535]$$

DICOM Protocol

- **Basic DICOM protocol actions**
 - **Test the connection** between two devices (C-Echo).
 - **Send images** from the local imaging device to a remote device (C-Store).
 - **Search the content** of a remote device (C-Find).
 - **Retrieve images** from a remote device (C-Move).
- SCU / SCP / Commands



Client (SCU)



dicominfo('file.dcm')

Filename: [1x55char]
FileModDate: '27-maig-200823: 25: 47'
FileSize: 17042728
Format: 'DICOM'
FormatVersion: 3
Width: 2560
Height: 3328
BitDepth: 12
ColorType: 'grayscale'
FileMetaInformationGroupLength: 212
FileMetaInformationVersion: [2x1uint8]
MediaStorageSOPClassUID: '1.2.840.10008.5.1.4.1.1.1.2'
MediaStorageSOPInstanceUID: [1x55char]
TransferSyntaxUID: '1.2.840.10008.1.2'
ImplementationClassUID: '1.3.6.1.4.1.5962.99.2'
ImplementationVersionName: 'PIXELMEDJAVA001'
SourceApplicationEntityTitle: 'RXMAMOPACS'
SpecificCharacterSet: 'ISO_IR100'
ImageType: 'DERIVED\PRIMARY\RIGHT'
SOPClassUID: '1.2.840.10008.5.1.4.1.1.1.2'
SOPInstanceUID: [1x55char]
StudyDate: '20070803'
SeriesDate: '20070803'
AcquisitionDate: '20070803'
ContentDate: '20070803'
StudyTime: '120122.375000'
SeriesTime: '120334.390000'
AcquisitionTime: '120310.812000'

ContentTime: '120330.281000'
AccessionNumber: '660471'
Modality: 'MG'
PresentationIntentType: 'FOR PRESENTATION'
Manufacturer: 'SIEMENS'
InstitutionName: 'HOSPITAL JOSEP TRUETA'
InstitutionAddress: [1x27char]
ReferringPhysicianName: [1x1struct]
StationName: 'YBFR015811'
SeriesDescription: 'Mamobilateral, Diagnose'
InstitutionalDepartmentName: 'Department'
PerformingPhysicianName: [1x2struct]
OperatorName: [1x1struct]
ManufacturerModelName: 'Mammomat Novation DR'
ReferencedPerformedProcedureStepSequence: [1x1struct]
DerivationDescription: [1x42char]
SourceImageSequence: [1x1struct]
AnatomicRegionSequence: [1x1struct]
PatientName: [1x1struct]
PatientID: '0000'
PatientBirthDate: '19560101'
PatientSex: 'F'
PatientAge: '051Y'
ContrastBolusAgent: ''
BodyPartExamined: 'BREAST'
KVP: 25
DeviceSerialNumber: '1036'
SoftwareVersion: [1x47char]

dicominfo('file.dcm')... cont

DistanceSourceToDetector: 650

DistanceSourceToPatient: 636

TableAngle: 0

TableType: 'NONE'

FieldOfViewShape: 'RECTANGLE'

FieldOfViewDimensions: [2x1double]

ExposureTime: 410

XrayTubeCurrent: 147

Exposure: 60

ExposureInuAs: 60099

RectificationType: 'CONST POTENTIAL'

ImagerPixelSpacing: [2x1double]

Grid: 'FOCUSEDPARALLEL'

FocalSpot: 0.3000

AnodeTargetMaterial: 'TUNGSTEN'

BodyPartThickness: 27

CompressionForce: 112

DateOfLastCalibration: '18590101'

TimeOfLastCalibration: '000000.000000'

AcquisitionDeviceProcessingDescription: ''

AcquisitionDeviceProcessingCode: ''

TomoLayerHeight: 0

PositionerType: 'MAMMOGRAPHIC'

PositionerPrimaryAngle: 0

PositionerSecondaryAngle: 0

DetectorPrimaryAngle: 0

DetectorSecondaryAngle: 0

ShutterShape: 'RECTANGULAR'

ShutterLeftVerticalEdge: 0

ShutterRightVerticalEdge: 2560

ShutterUpperHorizontalEdge: 0

ShutterLowerHorizontalEdge: 3328

ViewPosition: 'CC'

Sensitivity: 0

DetectorConditionsNominalFlag: 'YES'

DetectorTemperature: 30

DetectorType: 'DIRECT'

DetectorConfiguration: 'AREA'

DetectorDescription: ''

DetectorMode: 'Normalmode'

DetectorID: 'MP1775'

DateOfLastDetectorCalibration: '20070724'

TimeOfLastDetectorCalibration: '082552.265000'

ExposuresOnDetectorSinceLastCalibration: 23702

ExposuresOnDetectorSinceManufactured: 57839

DetectorTimeSinceLastExposure: 920.1090

DetectorActiveTime: 410

DetectorActivationOffsetFromExposure: 0

DetectorBinning: [2x1double]

DetectorElementPhysicalSize: [2x1double]

DetectorElementSpacing: [2x1double]

DetectorActiveShape: 'RECTANGLE'

DetectorActiveDimensions: [2x1double]

DetectorActiveOrigin: [2x1double]

FieldOfViewOrigin: [2x1 double]

FieldOfViewRotation: 0

FieldOfViewHorizontalFlip: 'NO'

FilterMaterial: 'RHODIUM'



dicominfo('file.dcm')... cont

FilterThicknessMinimum: 0.0500
 FilterThicknessMaximum: 0.0500
ExposureControlMode: 'AUTOMATIC'
 ExposureControlModeDescription: [1x55char]
 ExposureStatus: 'NORMAL'
 ExposureTimeInuS: 410000
 XrayTubeCurrentInuA: 1.4658e+005
StudyInstanceUID: [1x55char]
SeriesInstanceUID: [1x55char]
 StudyID: '1'
 SeriesNumber: 2
 AcquisitionNumber: 1
 InstanceNumber: 1
 PatientOrientation: 'P\L'
ImageLaterality: 'R'
 ImageComments: ''
 SamplesPerPixel: 1
 PhotometricInterpretation: 'MONOCHROME2'
 Rows: 3328
 Columns: 2560
 BitsAllocated: 16
 BitsStored: 12
 HighBit: 11
 PixelRepresentation: 0
 QualityControlImage: 'NO'
 BurnedInAnnotation: 'NO'
PixelIntensityRelationship: 'LOG'
 PixelIntensityRelationshipSign: 1



WindowCenter: [2x1double]
WindowWidth: [2x1double]
 RescaleIntercept: 0
 RescaleSlope: 1
 RescaleType: 'US'
 ImplantPresent: 'NO'
 PartialView: 'NO'
 WindowCenterWidthExplanation: 'WINDOW1\WINDOW2'
 LossyImageCompression: '00'
 Private_0029_10xx_Creator: 'SIEMENS MEDCOM HEADER'
 Private_0029_11xx_Creator: 'SIEMENS MEDCOM HEADER2'
Private_0029_1031: [12x1uint8]
Private_0029_1032: [4x1uint8]
Private_0029_1033: [4x1uint8]
Private_0029_1034: [12x1uint8]
 Private_0029_1160: [4x1uint8]
 StudyPriorityID: 'MED'
 PerformedProcedureStepStartDate: '20070803'
 PerformedProcedureStepStartTime: '120122.375000'
 PerformedProcedureStepID: 'MG20070803120122'
 EntranceDose: 0
 DistanceSourceToEntrance: 609
 OrganDose: 0.0070
OrganExposed: 'BREAST'
 AcquisitionContextSequence: [1x1struct]
 EntranceDoseInmGy: 2
 Private_0041_10xx_Creator: 'SIEMENSMEDSPDXMGWHAWS1'
 Private_0041_1002: [10x1uint8]

CalibrationImage: 'NO'
 ViewCodeSequence: [1x1struct]
 HistogramSequence: [1x1struct]
 Private_0095_10xx_Creator: 'SIENET'
 Private_0095_100c: [4x1uint8]
 PresentationLUTShape: 'IDENTITY'

`dicomread('file.dcm')`



dicomread('file.dcm')



To have in mind...

- Data protection! We are dealing with personal data!

http://plastimatch.org/dicom_comparison.html

- Image formats. 2D, 3D, 4D.
- Information visualization
- Bit depth: 8, 12, 16?
- Dimensionality, resolution,...
- Other formats: http://www.itk.org/Wiki/ITK/File_Formats
 - nrrd (Nearly Raw raster data <http://teem.sourceforge.net/nrrd/>)
 - mhd (meta image)
 - Analyze, etc...

DICOM Viewers

- [OsiriX](#)
- [Dicomworks](#)
- [MicroDicom](#)
- [JiveX DICOM Viewer](#)
- [3DimViewer](#)
- [Mango](#)
- [Escape](#)
- [IRFANVIEW](#)
- [RadiAnt](#)
- [HOROS](#)
- [Itk Snap](#)

To know more

- <https://www.dicomstandard.org/> Official Site
- <https://en.wikipedia.org/wiki/DICOM>
- [Itk.org](http://itk.org)
- Insight into images